



Exclusive distributors of Model-Med OB/GYN simulators in North America  
Visit us online at [parmedsim.com](http://parmedsim.com) or find us at the following conferences

### Conference Schedule

IMSH 2022 Conference

January 15 to 19 – Los Angeles, CA

[imsh2022.org](http://imsh2022.org)

ACOGG Annual Conference

April 3 to 8 – San Antonio, TX

[acog.org/events](http://acog.org/events)

SMFM Annual Pregnancy Meeting

Jan 31 to Feb 5 – Orlando, FL

[smfm.org/2022](http://smfm.org/2022)

ACOG Annual Meeting

May 6 to 8 – San Diego, CA

[annualmeeting.acog.org](http://annualmeeting.acog.org)

CREOG & APGO Annual Meeting

March 9 to 12 – Orlando, FL

[apgo.org/2022caam](http://apgo.org/2022caam)

ACNM Annual Meeting

May 22 to 26 – Chicago, IL

[annualmeeting.midwife.org](http://annualmeeting.midwife.org)

Paradigm Medical Systems' full line of Model-med trainers are constructed with a proprietary polymer and designed with an emphasis on accurate internal anatomy to bring incredible tactile realism to your simulations. Remarkable polymer elasticity and durability coupled with the rigidity of accurate internal skeletal structures creates an unparalleled experience, enabling lifelike application of maneuvers and realistic outcomes intrinsic to the model's design.

Model-med trainers do not require electricity or rely on electronic components, reducing complexity and increasing reliability. By focusing on innovative design and materials, Sophie and Her Mum is able to closely mimic maternal and fetal anatomy both in terms of corporal consistency and structural rigidity. The same polymer is molded over the maternal model's internal bony pelvis and the fetal model's fully articulated skeletal system, achieving a realistically snug birth canal and a full range of fetal motion. The accurate maternal anatomy and articulated fetal model enables mimicry of any feasible presentation.

#### **Paradigm Medical Systems**

North American Distributor

Web: [parmedsim.com](http://parmedsim.com)

Phone: (971) 570-4244

Email: [bridget@parmedsim.com](mailto:bridget@parmedsim.com)

#### **Model-med** (contact if outside the Americas)

Manufacturer

Web: [modelmed.com.au](http://modelmed.com.au)

Phone: +61 438 353 910

Email: [sales@modelmed.com.au](mailto:sales@modelmed.com.au)